EMCal LEDs Testing


Apurva Narde

University of Illinois Urbana-Champaign


April 7, 2023

Introduction

- Each block has four towers (readout channels) with four SiPMs (Silicon Photomultipliers) each, resulting in a total of:
 - ▶ 24,576 readout channels
 - ▶ 98,304 SiPMs
 - ▶ 16.5 million fibers.
- ➤ There are 1,560 scintillating fibers extending along the length direction are embedded in the block.
- ► Each tower represents 1 ADC (Analog-to-digital converter) readout channel and is equipped with a light guide coupled to 4 SiPMs that collect the light from the fibers.




SiPM PCB daughterboard 4 SiPMs per tower



4 towers per block with channel ordering

Readout Channels and Interface Board

- ▶ There are $96 \times 4 = 384$ distinctive ADC channels in each sector.
- No physical boundary between neighboring readout towers within a block.
- Channels are organized into groups of 64, each group (6 total) controlled by an interface board (IB).



LEDs Testing


- ightharpoonup Electronics to readout chain: Signal ightharpoonup SiPM ightharpoonup Readout.
- ► LED's go in place of "Signal".
- ► Test pulse circumvents the SiPM and goes directly to readout.
- ► Test pulses act as the ideal response whereas LEDs are used to measure the performance of the SiPMs.

Test Pulse Data

- ► Location:
 - /sphenix/lustre 01/sphnx pro/commissioning/emcal/calib/calib-00005027-0000.prdf
- Events: 10,999.
 - Each event contains data for time, adc, and pedestal for each channel.
- Currently data is only for 4 sectors.
- ▶ Total Channels: 1,536 = 4 sectors \times 384 channels per sector.
 - ▶ Note: Event 0 has data for 1,152 channels and event 1 has data for 1,408 channels.
 - Remaining events have data for 1,536 channels.
- Packets: 6009 to 6016 (8 packets total).
- ightharpoonup 1 sector = 2 packets, 1 packet = 3 IBs, 1 IB = 64 channels.


Time







 Adc






Ped Pedestal 10⁷ Counts 10⁶ 10⁵ 10⁴ 10^{3} 10^2 10 1900 Ped 1300 1400 1500 1600 1700 1800




Time vs Channel






Adc vs Channel






Pedestal vs Channel



Adc vs Time



Pedestal vs Time

